

# OUR PRODUCTS



## EXAMINATION MICROSCOPES

Magnify skin and hair from 1 to 500 times while displaying the image on a TV, computer screen, iPad, Android tablet, LCD panel or projector so you and the client can see. The microscopes are very easy to use. Simply touch the area of the skin, scalp, or hair to be viewed and see the image on the screen. Skin, scalp, and hair conditions are easier to see at various magnifications and dramatically show the cause of the condition and further motivate the client to use your products and services to solve the problem.




## LENSES

All of our microscopes offer several lenses or focal tips and a light source to greatly enhance the image of the object you are viewing. Although most are standard lenses, polarized tips and lenses are available along with lenses with internal lights. Using the various lenses and tips allow for the widest magnification possible for each microscope and is suitable for most settings.



## SKIN AND HAIR KITS

Kits contain all the equipment and materials needed for a complete examination of a client's skin or hair, consequently raising their awareness of their condition and their desire to correct the problem with your products and/or services.



## EXAMINATION AND ANALYSIS SOFTWARE

CompareView and ComparView Hair is an easy to use but powerful software program. It assists users with multiple on-screen image comparisons, client management, basic image editing, and image analysis for alopecia, hair density, and hair width following the Norman scale. CompareView software is compatible with Windows computers at this time.